

Product Name: Recombinant Human YTHDF1 Protein
Catalog #: PEH10638



Summary

Name	YTHD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BYJ9
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDTLSKAPGM NSLEQGMVGL KIGDVSSAV KTVGSVSSV ALTGVLSGNG GTNVNMPVSK PTSWAAIASK PAKPQPKMKT KSGPVMGGGL PPPPIKHNMD IGTWDNKGPV PKAPVPQQAP SPQAAPQPPQ VAQPLPAQPP ALAQPQYQSP QQPPQTRWVA PRNRNAAFQ SGGAGSDSNS PGNVQPNSAP SVESHPVLEK LKAHSYNPK EFEWNLKSGR VFIKSYSED DIHRSIKYSI WCSTEHGNKR LDSAFRCMSS KGPVYLLFSV NGSFHFCGVA EMKSPVDYGT SAGVWSQDKW KGFQDVQWIF VKDVPNNQLR HIRLENNDNK PVTNSRDTQE VPLEKAKQVL KIISSYKHTT SIFDDFAHYE KRQEEEEVVR

image



Background

Alternative Names YTH domain-containing family protein 1 (Dermatomyositis associated with

Product Name: Recombinant Human YTHDF1 Protein
Catalog #: PEH10638



cancer putative autoantigen 1) (DACA-1)

Note

For Research Use Only , Not for Diagnostic Use