


Product Name: Recombinant Human RPL8 Protein
Catalog #: PEH10640



Summary

Name	RL8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P62917
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKGIVKDIIH DPGRGAPLAK VVFRDPYRFK KRTELFIAAE GIHTGQFVYC GKKAQLNIGN VLPVGTMPPEG TIVCCLEEKP GDRGKLARAS GNYATVISHN

image



Background

Alternative Names 60S ribosomal protein L8 (Large ribosomal subunit protein uL2)

Note

For Research Use Only , Not for Diagnostic Use