


Product Name: Recombinant Human CYTH1 Protein
Catalog #: PEH10642



Summary

Name	CYH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15438
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLTAEERQEL ENIRRRKQEL LADIQLKDE IAEVANEIEN LGSTEERKNM QRNKQVAMGR KKFNMDPKKG IQFLIENDLL KNTCEDIAQF LYKGEGLNKT AIGDYLGERD EFNIQVLHAF VELHEFTDLN LVQALRQFLW SFRLPGEAQK IDRMMEAFQA RYCQCNGGVF QSTDTCYVLS FAIIMLNTSL HNPNVKDKPT VERFIAMNRG INDGGDLPEE LLRNLYESIK NEPFKIPEDD GNDLTHTFN PDREGWLLKL GGRVKTWKR RWFILTDNCL YFEYTTDKE PRGIIPLENL SIREVEDSKK PNCFELYIPD NKDQVIKACK TEADGRVVEG NHTVYRISAP TPEEKEEWIK CIKAAISRDP

image



Background

Alternative Names Cytohesin-1 (PH, SEC7 and coiled-coil domain-containing protein 1) (SEC7)

Product Name: Recombinant Human CYTH1 Protein
Catalog #: PEH10642



homolog B2-1)

Note

For Research Use Only , Not for Diagnostic Use