


Product Name: Recombinant Human BEX2 Protein
Catalog #: PEH10644



Summary

Name	BEX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXY8
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MESKEERALN NLIVENVNQE NDEKDEKEQV ANKGEPLALP LNVSEYCVPR GNRRRFRVRQ PILQYRWDIM HRLGEPQARM REENMERIGE EVRQLMEKLR EKQLSHSLRA VSTDPHHDH HDEFCLMP

image



Background

Alternative Names Protein BEX2 (Brain-expressed X-linked protein 2) (hBex2)

Note

For Research Use Only , Not for Diagnostic Use