


## Summary

<b>Name</b>	SPCS_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9HD40
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	21 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DQKSCFKSMI TAGFEPVVIE NVLEGDELRT DLKAVEAKVQ ELGPDCILCI HSTTSCFAPR VPDRLEELAV ICANYDIPHI VNNAYGVQSS KCMHLIQQGA RVGRIDAFVQ SLDKNFMVPV GGAIAGFND SFIQEISKMY PGRASASPSL DVLITLLSLG SNGYKLLKE

## image



## Background

<b>Alternative Names</b>	O-phosphoseryl-tRNA(Sec) selenium transferase (EC 2.9.1.2) (Liver-pancreas antigen) (LP) (SLA-p35) (SLA/LP autoantigen) (Selenocysteine synthase) (Sec synthase) (Selenocysteinyl-tRNA(Sec) synthase) (Sep-tRNA:Sec-tRNA synthase) (SepSecS) (Soluble liver antigen) (SLA) (UGA suppressor tRNA-associated protein) (tRNA(Ser/Sec)-associated antigenic protein)
--------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Product Name: Recombinant Human SEPSECS Protein**  
**Catalog #: PEH10645**

---



**Note**

For Research Use Only , Not for Diagnostic Use