


Product Name: Recombinant Human SEPSECS Protein
Catalog #: PEH10646



Summary

Name	SPCS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HD40
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TAGFEPVVIE NVLEGDELRT DLKAVEAKVQ ELGPDCLCI HSTTSCFAPR VPDRLEELAV ICANYDIPHI VNNAYGVQSS KCMHLIQQGA RVGRIDAFVQ SLDKNFMVPV GGAIAGFND SFIQEISKMY PGRASAPSL DVLITLLSLG SNGYKLLKE RKEMFSYLSN QIKKLSEAYN

image



Background

Alternative Names	O-phosphoseryl-tRNA(Sec) selenium transferase (EC 2.9.1.2) (Liver-pancreas antigen) (LP) (SLA-p35) (SLA/LP autoantigen) (Selenocysteine synthase) (Sec synthase) (Selenocysteinyl-tRNA(Sec) synthase) (Sep-tRNA:Sec-tRNA synthase) (SepSecS) (Soluble liver antigen) (SLA) (UGA suppressor tRNA-associated protein) (tRNA(Ser/Sec)-associated antigenic protein)
--------------------------	--

Product Name: Recombinant Human SEPSECS Protein
Catalog #: PEH10646



Note

For Research Use Only , Not for Diagnostic Use