


Summary

Name	TRIM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9C040
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDKDGELCKT GNAYLTAELS TPDGSVADGE ILDNKNGTYE FLYTVQKEGD FTLSRLYDQ HIRGSPFKLK VIRSADVST TEGVKRRVKS PGSGHVKQKA VKRPASMYST GKRKENPIED DLIFRVGTKG

image



Background

Alternative Names	Tripartite motif-containing protein 2 (EC 2.3.2.27) (E3 ubiquitin-protein ligase TRIM2) (RING finger protein 86) (RING-type E3 ubiquitin transferase TRIM2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRIM2 Protein
Catalog #: PEH10656

