


Product Name: Recombinant Human NUP85 Protein
Catalog #: PEH10668



Summary

Name	NUP85_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BW27
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALPLLEQKQV IFSAEQTYEL MRCLEDLTSR RPVHGESDTE QLQDDDIETT KVEMLRSLA RNLARAIRE

image



Background

Alternative Names Nuclear pore complex protein Nup85 (85 kDa nucleoporin) (FROUNT)
(Nucleoporin Nup75) (Nucleoporin Nup85) (Pericentrin-1)

Note

For Research Use Only , Not for Diagnostic Use