


Product Name: Recombinant Human PARVA Protein
Catalog #: PEH10669



Summary

Name	PARVA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NVD7
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPKSPTPKSP PSRKKDDSL GKLGGTLARR KKAKEVSELQ EEGMNAINLP LSPIPFELDP EDTMLEENEV RTMVDPNSRS

image



Background

Alternative Names	Alpha-parvin (Actopaxin) (CH-ILKBP) (Calponin-like integrin-linked kinase-binding protein) (Matrix-remodeling-associated protein 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use