


Product Name: Recombinant Human DYNLRB1 Protein
Catalog #: PEH10673



Summary

| | |
|------------------------------|---|
| Name | DLRB1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9NP97 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 12 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | MAEVEETLKR LQSQKGVQGI IVVNTGIIPI KSTMDNPTTT QYASLMHSFI LKARSTVRDI DPQNDLTLFLR IRSKKNEIMV APDKDYFLIV IQNPTE |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Dynein light chain roadblock-type 1 (Bithoraxoid-like protein) (BLP) (Dynein light chain 2A, cytoplasmic) (Dynein-associated protein Km23) (Roadblock domain-containing protein 1) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use