

Product Name: Recombinant Human KATNA1 Protein
Catalog #: PEH10677



Summary

Name	KTNA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75449
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKGREEKNKS PAAVTEPETN KFDSTGYDKD LVEALERDII SQNPVVRWDD IADLVEAKKL LKEAVVLP MW

image



Background

Alternative Names Katanin p60 ATPase-containing subunit A1 (Katanin p60 subunit A1) (EC 5.6.1.1) (p60 katanin)

Note

For Research Use Only , Not for Diagnostic Use