


Product Name: Recombinant Human STAP1 Protein
Catalog #: PEH10682



Summary

Name	STAP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9ULZ2
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLREKKRRI ETEQSTSVEK EKEPTEDYVD VLNPMPACFY TVSRKEATEM LQKNPSLGNM ILRPGSDSRN YSITIRQEID IPRIKHYKVM SVGQNYTIEL EKPVTLPNLF SVIDYFVKET RGNLRPFICS TDENTGQEPS MEGRSEKLLK

image



Background

Alternative Names Signal-transducing adaptor protein 1 (STAP-1) (BCR downstream-signaling protein 1) (Docking protein BRDG1) (Stem cell adaptor protein 1)

Note

For Research Use Only , Not for Diagnostic Use