


Product Name: Recombinant Human PLCG1 Protein
Catalog #: PEH10691



Summary

Name	PLCG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P19174
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YGALYEGRNP GFYVEANPMP TFKCAVKALF DYKAQREDEL TFIKSAIIQN VEKQEGGWWR GDYGGKKQLW FPSNYVEEMV NPVALEPERE HLDENSPLGD LLRGVLDVPA CQIAIRPEGK NNRLFVFSIS MASVAHWSLD VAADSQEELQ DWVKKIREVA QTADARLTEG KIMERRKKIA

image



Background

Alternative Names	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1 (EC 3.1.4.11) (PLC-148) (Phosphoinositide phospholipase C-gamma-1) (Phospholipase C-II) (PLC-II) (Phospholipase C-gamma-1) (PLC-gamma-1)
--------------------------	--

Note

Product Name: Recombinant Human PLCG1 Protein
Catalog #: PEH10691



For Research Use Only , Not for Diagnostic Use