


Product Name: Recombinant Human RPS11 Protein
Catalog #: PEH10693



Summary

Name	RS11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P62280
Host	E.coli
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MADIQTERAY QKQPTIFQNK KRVLLGETGK EKLPRYYKNI GLGFKTPKEA IEGTYIDKKC PFTGNVSIRG RILSGVVTM KMQRITVIRR DYLYIRKYN RFEKRHNMS VHLSPCFRDV QIGDIVTGE CRPLSKTVRF NVLKVTKAAG TKKQFQKF

image



Background

Alternative Names 40S ribosomal protein S11 (Small ribosomal subunit protein uS17)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RPS11 Protein
Catalog #: PEH10693

