



Product Name: Recombinant Human LEF1 Protein
Catalog #: PEH10697



Summary

Name	LEF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJU2
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEIIPASNGH EVARQAQTSQ EPYHDKAREH PDDGKHPDGG LYNKGPSYSS YSGYIMMPNM NNDPYMSNGS LSPPIRPTSN KVPVVQPSHA VHPLTPLITY SDEHFSPGSH PSHIPSDVNS KQGMSRHPPA PDIFTFYPLS PGGVGQITPP LGWQGQPVYP ITGGFRQPYP SLSVDTSMS RFSHHMIPGP PGPHTTGIPH PAIVTPQVKQ EHPHTDSDLM HVKPQHEQRK EQEPKRPHIK

image



Background

Alternative Names	Lymphoid enhancer-binding factor 1 (LEF-1) (T cell-specific transcription factor 1-alpha) (TCF1-alpha)
--------------------------	--

Note

Product Name: Recombinant Human LEF1 Protein
Catalog #: PEH10697



For Research Use Only , Not for Diagnostic Use