


Product Name: Recombinant Human TTK Protein
Catalog #: PEH10699



Summary

Name	TTK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P33981
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QNSHFKEPLV SDEKSELII TDSITLKNKT ESSLLAKLEE TKEYQEPEVP ESNQKQWQSK RKSECINQNP AASSNHWQIP ELARKVNTEQ KHTTFEQPVF SVSKQSPPIIS

image



Background

Alternative Names	Dual specificity protein kinase TTK (EC 2.7.12.1) (Phosphotyrosine picked threonine-protein kinase) (PYT)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use