


Product Name: Recombinant Human INA Protein
Catalog #: PEH10704



Summary

Name	AINX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16352
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASDGLDLSQA AARTNEYKII RTNEKEQLQG LNDRFVAFIE KVHQLTQNR ALEAELAALR QRHAEPSRVG ELFQRELRDL RAQLEEASSA RSQALLERDG LAEEVQRLRA RCEEESRGRE GAERALKAAQ RDVDGATLAR LDLEKKVESL LDELAFVRQV HDEEVAELLA TLQASSQAAA EVDVTVAKPD LTSALREIRA QYESLAAKNL QSAEEWYKSK FANLNEQAAR STEAIRASRE EIHEYRRQLQ ARTIEIEGLR GANESLERQI LELEERHSAE VAGYQDSIGQ LENDLRNTKS EMARHLREYQ DLLNVKMALD IEIAAYRKLL

image



Background

Alternative Names	Alpha-interneixin (Alpha-Inx) (66 kDa neurofilament protein) (NF-66) (Neurofilament-66) (Neurofilament 5)
--------------------------	---

Product Name: Recombinant Human INA Protein
Catalog #: PEH10704



Note

For Research Use Only , Not for Diagnostic Use