


Product Name: Recombinant Human UBE2C Protein
Catalog #: PEH10707



Summary

Name	UBE2C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00762
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASQNRDPAA TSVAARKGA EPSGGAARGP VGKRLQQELM TLMMSGDKGI SAFPESDNLF KWVGTIHGAA GTVYEDLRYK LSLEFSGYP YNAPTVKFLT PCYHPNVDQT GNICLDILKE KWSALYDVRT ILLSIQSLG EPNIDSPLNT HAAELWKNPT AFKKYLQETY SKQVTSQEP

image



Background

Alternative Names	Ubiquitin-conjugating enzyme E2 C (EC 2.3.2.23) ((E3-independent) E2 ubiquitin-conjugating enzyme C) (EC 2.3.2.24) (E2 ubiquitin-conjugating enzyme C) (UbcH10) (Ubiquitin carrier protein C) (Ubiquitin-protein ligase C)
--------------------------	--

Note

Product Name: Recombinant Human UBE2C Protein
Catalog #: PEH10707



For Research Use Only , Not for Diagnostic Use