


Product Name: Recombinant Human CDKN2C Protein
Catalog #: PEH10710



Summary

Name	CDN2C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42773
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAEPWGNELA SAAARGDLEQ LTSLLQNNVN VNAQNGFGRT ALQVMKLGNP EIARRLLLRG ANPDLKDRTG FAVIHDAARA GFLDTLQTLL EFQADVNIED NEGNLPLHLA AKEGHLRVVE FLVKHTASNV GHRNHKGDTA CDLARLYGRN EVSLSMQANG AGGATNLQ

image



Background

Alternative Names	Cyclin-dependent kinase 4 inhibitor C (Cyclin-dependent kinase 6 inhibitor) (p18-INK4c) (p18-INK6)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CDKN2C Protein
Catalog #: PEH10710

