


Product Name: Recombinant Human GNS Protein
Catalog #: PEH10712



Summary

Name	GNS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P15586
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRWQTLLSVD DLVEKLVKRL EFTGELNNTY IFYTSNNGYH TGQFSLPIDK RQLYEFDIKV PLLVRGPGIK PNQTSKMLVA NIDLGPITLD IAGYDLNKTK MDGMSLLPIL RGASNLTWRS DVLVEYQGEG

image



Background

Alternative Names N-acetylglucosamine-6-sulfatase (EC 3.1.6.14) (Glucosamine-6-sulfatase) (G6S)

Note

For Research Use Only , Not for Diagnostic Use