


Product Name: Recombinant Human SENP8 Protein
Catalog #: PEH10717



Summary

Name	SENP8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96LD8
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDPVVLSYMD SLLRQSDVSL LDPPSWLNDH IIGFAFEYFA NSQFHDCSDH VSFISPEVTQ FIKCTSNPAE IAMFLEPLDL PNKRVVFLAI NDNSNQAAGG THWSLLVYLQ DKNSFFHYDS HRSNSVHAK QVAEKLEAFL GRKGDKLAFLV EEKAPAQQNS YDCGMYVICN TEALCQNFFR QQTESLLQLL TPAYITKKRG EWKDLITTLA KK

image



Background

Alternative Names Sentrin-specific protease 8 (EC 3.4.22.-) (Deneddylase-1) (NEDD8-specific protease 1) (Protease, cysteine 2) (Sentrin/SUMO-specific protease SENP8)

Note

Product Name: Recombinant Human SENP8 Protein
Catalog #: PEH10717



For Research Use Only , Not for Diagnostic Use