


Product Name: Recombinant Human UBE2M Protein
Catalog #: PEH10720



Summary

Name	UBC12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P61081
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKEESAGGT KGSSKKASAA QLRIQKDINE LNLPKTCDIS FSDPDDLNF KLVICPDEGF YKSGKFVFSF

image



Background

Alternative Names	NEDD8-conjugating enzyme Ubc12 (EC 2.3.2.34) (NEDD8 carrier protein) (Ubiquitin-conjugating enzyme E2 M)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use