


**Product Name: Recombinant Human DNAJB9 Protein**  
**Catalog #: PEH10722**



## Summary

<b>Name</b>	DNJB9_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UBS3
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	18 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LGVPKSASER QIKKAFHKLA MKYHPDKNKS PDAEAKFREI AEAYETLSDA NRRKEYDTLG HSAFTSGKGQ RGS GSSFEQS FNFNFDLDFK DFGFFGQNQN TGSKKRFENH FQTRQDGGSS RQRHHFQEFS FGGGLFDDMF

## image



## Background

<b>Alternative Names</b>	DnaJ homolog subfamily B member 9 (Endoplasmic reticulum DNA J domain-containing protein 4) (ER-resident protein ERdj4) (ERdj4) (Microvascular endothelial differentiation gene 1 protein) (Mdg-1)
--------------------------	--

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human DNAJB9 Protein**  
**Catalog #: PEH10722**

---

