


Summary

Name	AADAT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5Z0
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AYSGTLQSLH PLGCNIINVA SDESGIVPDS LRDILSRWKP EDAKNPQKNT PKFLYTPNG NNPTGNSLTS

image



Background

Alternative Names	Kynurenine/alpha-aminoadipate aminotransferase, mitochondrial (KAT/AadAT) (2-aminoadipate aminotransferase) (2-aminoadipate transaminase) (EC 2.6.1.39) (Alpha-aminoadipate aminotransferase) (AadAT) (Kynurenine aminotransferase II) (Kynurenine--oxoglutarate aminotransferase II) (Kynurenine--oxoglutarate transaminase 2) (EC 2.6.1.7) (Kynurenine--oxoglutarate transaminase II)
--------------------------	---

Product Name: Recombinant Human AADAT Protein
Catalog #: PEH10723



Note

For Research Use Only , Not for Diagnostic Use