


Product Name: Recombinant Human PPP1R3C Protein
Catalog #: PEH10727



Summary

Name	PPR3C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UQK1
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPYDEFQRRH FVNKLKPLKS CLNIHKAKS QNDWKCSHNQ AKKRVPFADS KGLSLTAIHV FSDLPEEPAW DLQFDLLDLN DISSALKHHE EKNLILDFPQ

image



Background

Alternative Names	Protein phosphatase 1 regulatory subunit 3C (Protein phosphatase 1 regulatory subunit 5) (PP1 subunit R5) (Protein targeting to glycogen) (PTG)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use