


Product Name: Recombinant Human LGALS2 Protein
Catalog #: PEH10731



Summary

Name	LEG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05162
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTGELEVKNM DMKPGSTLKI TGSADGTDG FVINLGQGTG KLNLFNPRF SESTIVCNSL DGSNWGQEQR EDHLCFSPGS EVKFTVTFES DKFKVKLPDG HELTFPNRLG HSHLSYLSVR GGFNMSSFKL KE

image



Background

Alternative Names	Galectin-2 (Gal-2) (Beta-galactoside-binding lectin L-14-II) (HL14) (Lactose-binding lectin 2) (S-Lac lectin 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use