


Product Name: Recombinant Human FKBP2 Protein
Catalog #: PEH10735



Summary

Name	FKBP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P26885
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATGAEGKRKL QIGVKKRVDH CPIKSRKGDV LHMHYTGKLE DGTEFDSSLP QNQPFVFSLG TGQVIKGWDQ GLLGMCEGEK RKLVIPSELG YGERGAPPKI PGGATLVFEV ELLKIERRTE L

image



Background

Alternative Names	Peptidyl-prolyl cis-trans isomerase FKBP2 (PPIase FKBP2) (EC 5.2.1.8) (13 kDa FK506-binding protein) (13 kDa FKBP) (FKBP-13) (FK506-binding protein 2) (FKBP-2) (Immunophilin FKBP13) (Rotamase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FKBP2 Protein
Catalog #: PEH10735

