


Product Name: Recombinant Human NCOR2 Protein
Catalog #: PEH10736




Summary

Name	NCOR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y618
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPSRLDRGRE DSLPKGHVIY EGKKGHVLSY EGGMSVTQCS KEDGRSSSGP PHETAAPKRT YDMMEGRVGR AISSASIEGL MGRAIPPERH SPHHLKEQHH IRGSITQGIP RSYVEAQEDY LRREAKLLKR EGTPPPPPPS RDLTEAYKTQ ALGPLKLKPA HEGLVATVKE AGRSIHEIPR EELRHTPELP LAPRPLKEGS ITQGTPLKYD

image



SDS



WB

Background

Alternative Names	Nuclear receptor corepressor 2 (N-CoR2) (CTG repeat protein 26) (SMAP270) (Silencing mediator of retinoic acid and thyroid hormone receptor) (SMRT) (T3 receptor-associating factor) (TRAC) (Thyroid-, retinoic-acid-receptor-associated corepressor)
--------------------------	---

Product Name: Recombinant Human NCOR2 Protein
Catalog #: PEH10736



Note

For Research Use Only , Not for Diagnostic Use