


Product Name: Recombinant Human POGLUT1 Protein
Catalog #: PEH10742



Summary

Name	PGLT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8NBL1
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KWKVFIDQIN RSLNYPECS SQNCSCYHGV IEEDLTPFRG GISRKMMMAEV VRRKLGTHYQ ITKNRLYREN DCMFPSRCSG VEHFIEVIG RLPDMMEMVIN VRDYPQVPKW MEPAIPVFSF SKTSEYHDIM

image



Background

Alternative Names	Protein O-glucosyltransferase 1 (EC 2.4.1.376) (CAP10-like 46 kDa protein) (hCLP46) (KTEL motif-containing protein 1) (Myelodysplastic syndromes relative protein) (O-glucosyltransferase Rumi homolog) (hRumi) (Protein O-xylosyltransferase POGLUT1) (EC 2.4.2.63)
--------------------------	--

Note

Product Name: Recombinant Human POGlut1 Protein
Catalog #: PEH10742



For Research Use Only , Not for Diagnostic Use