


Product Name: Recombinant Human MAP3K21 Protein
Catalog #: PEH10745



Summary

Name	M3K21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TCX8
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QESFHSMQDD WKLEIQMFD ELRTKEKELR SREEELTRAA LQKKSQEELL KRREQQLAER EIDVLERELN ILIFQLNQEY PKVKKRKGKF KRSRLKLKDG HRISLPSEFQ HKITVQASPN

image



Background

Alternative Names	Mitogen-activated protein kinase kinase kinase 21 (EC 2.7.11.25) (Mitogen-activated protein kinase kinase kinase MLK4) (Mixed lineage kinase 4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use