


Product Name: Recombinant Human MAGEB1 Protein
Catalog #: PEH10750



Summary

Name	MAGB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P43366
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AREKRRKARE ETQGLKVAHA TAAEKEECPS SSPVLGDTPT SSPAAGIPQK PQGAPPTTTA AAVSCTESD EGA KCQGEEN ASFSQATTST

image



Background

Alternative Names	Melanoma-associated antigen B1 (Cancer/testis antigen 3.1) (CT3.1) (DSS-AHC critical interval MAGE superfamily 10) (DAM10) (MAGE-B1 antigen) (MAGE-XP antigen)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use