


Product Name: Recombinant Human CCDC57 Protein
Catalog #: PEH10753



Summary

Name	CCD57_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2TAC2
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WRALQAHRTQ LQEALQDTR SQLEEAQGKL RCLQEDFVYN LQVLEERDLE LERYDAAFAQ AREWEEARRA EVSELKIEAA KLRQALAREA RKVEELQQQQ QLAFQEHRLE LERVHSDKNG EIDHHREQYE NLKWTLERKL EELDGEALQ RQELLLEFES KMRKREHEFR LQADNMSNTA LSRELKVKLL HKELEALKEA GAKAAESLQR AEATNAELER KLQSRAGELQ DLEAMSRARV KDLEDKLHSV QLTRKKEEET FKRKHEELDR LAREKDAVLV AVKGAHVEQL QELQTRVLEL QAHCETLEAQ LRRAEWRQAD TAKEKDAAID QLREDASTVK SAWDAQIAQL SKEMVSRDLQ IQTLQEEEVK LKAQVARSQQ DIERYKQQLS LAVERERSLE RDQVQLGLDW QRRCDIERD QIQKSEALIQ GLSMAKSQVA AKLQETEQAL QEQEVVLKAV

image



Background

Product Name: Recombinant Human CCDC57 Protein
Catalog #: PEH10753



Alternative Names

Coiled-coil domain-containing protein 57

Note

For Research Use Only , Not for Diagnostic Use