


Product Name: Recombinant Human CMC4 Protein
Catalog #: PEH10755



Summary

Name	CMC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56277
Host	E.coli
Species	Human
Molecular Mass	9 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPQKDPCQKQ ACEIQKCLQA NSYMESKCQA VIQELRKCCA QYPKGRSVVC SGFEKEEEEN LTRKSASK

image



Background

Alternative Names	Cx9C motif-containing protein 4 (Mature T-cell proliferation 1 neighbor protein) (Mature T-cell proliferation-1 type A) (MTCP-1 type A) (Protein p8 MTCP-1) (p8MTCP1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use