


Product Name: Recombinant Human MRPS11 Protein
Catalog #: PEH10760



Summary

Name	RT11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P82912
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTICTGARQL QDAAAKQKVE QNAAPSHTKF SIYPPIPGEE SSLRWAGKKF EEIPIAHKA SHNNTQIQVV SASNEPLAFA SCGTEGFRNA KKG TGIAAQT AGIAAAARAK QKGVIHIRVV VKGLGPGRSL AMHGLIMGGL EVISITDNTF

image



Background

Alternative Names	28S ribosomal protein S11, mitochondrial (MRP-S11) (S11mt) (Cervical cancer proto-oncogene 2 protein) (HCC-2) (Mitochondrial small ribosomal subunit protein uS11m)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MRPS11 Protein
Catalog #: PEH10760

