


Product Name: Recombinant Human CCDC180 Protein
Catalog #: PEH10764



Summary

Name	CC180_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P1Z9
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NDQEMESFIT EEVLGQQKKS PLHAKMDESK EGSIQGLEEM QVEREGSLNP SLNEENVKGQ GEKKEESEE DEKEEEEEEE KLEEEKEEKE AQEEQESLSV GEEEDKEEGL EEIYYEDMES FTISSGNTYF VFVPLEEEHC RKSHSTFSAM FINDTSSAKF IEQVTIPSRL ILEIKKQLRA GFFEHLKWF DQCSLNTRVT VATKINELDS ELELHLHLHQ PRAQQIEKDI HNVRAAELL HQEQLDSHCA GVTETLKKKR LMFCQFQEEQ NVRSKNFRLK IYDMEHIFLN ATRSQKLVTL SNTLHQELLS YVDVTQVSLR SFRQYLEESL GKLYRSNIEF IKHCRLFSEG GNFSPKEINS LCSRLEKEAA RIELVESVIM LNMEKLENEY LDQANDVINK FESKFHNLSV DLIFIEKIQR LLTNLQVKIK

image



Background

Product Name: Recombinant Human CCDC180 Protein
Catalog #: PEH10764



Alternative Names

Coiled-coil domain-containing protein 180

Note

For Research Use Only , Not for Diagnostic Use