


Product Name: Recombinant Human MRPL43 Protein
Catalog #: PEH10765



Summary

Name	RM43_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N983
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RYVQQLQRLS FVSRDGASS RGAREFVERE VIDFARRNPG VVIYVNSRPC CVPRVVAEYL NGAVREESIH CKSVEEISTL VQKLADQSGL DVIRIRKPFH TDNPSIQGQW HPFTNKPTTF RGLRPREVQD PAPAQDTGLR LSAVAPQILL PGWPDPPDLP

image



Background

Alternative Names	39S ribosomal protein L43, mitochondrial (L43mt) (MRP-L43) (Mitochondrial large ribosomal subunit protein mL43) (Mitochondrial ribosomal protein bMRP36a)
--------------------------	---

Note

Product Name: Recombinant Human MRPL43 Protein
Catalog #: PEH10765



For Research Use Only , Not for Diagnostic Use