


Product Name: Recombinant Human BCL2L2 Protein
Catalog #: PEH10766



Summary

Name	B2CL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92843
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RALVADFVGY KLRQKGYVCG AGPGEGPAAD PLHQAMRAAG DEFETFRRT FSDLAAQLHV TPGSAQQRFT QVSDELFQGG PNWGRVLAFF VFGAALCAES VNKEMEPLVG QVQEWVAYL ETQLADWIHS SGGWAEFTAL YGDGALEEAR RLREGNWASV

image



Background

Alternative Names Bcl-2-like protein 2 (Bcl2-L-2) (Apoptosis regulator Bcl-W)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BCL2L2 Protein
Catalog #: PEH10766

