


Product Name: Recombinant Human CA5B Protein
Catalog #: PEH10767



Summary

Name	CAH5B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y2D0
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CSLYTCTYKT RNRALHPLWE SVDLVPGGDR QSPINIRWRD SVYDPGLKPL TISYDPATCL HVWNNGYSL VEFEDSTDKS VIKGGPLEHN YRLQFHFHW GAIDAWGSEH TVDSKCFPAE LHLVHWNVR FENFEDAALE ENGLAVIGVF LKLGHKHKEL QKLVDLPSI KHKDALVEFG SFDPSCLMPT CPDYWTYSGS LTPPLSESV TWIHKQPVE VDHDQLEQFR TLLFTSEGEK EKRMVDNFRP LQPLMNRTVR SSFRHDYVLN VQAKPKPATS QATP

image



Background

Alternative Names	Carbonic anhydrase 5B, mitochondrial (EC 4.2.1.1) (Carbonate dehydratase VB) (Carbonic anhydrase VB) (CA-VB)
--------------------------	--

Product Name: Recombinant Human CA5B Protein
Catalog #: PEH10767



Note

For Research Use Only , Not for Diagnostic Use