


Product Name: Recombinant Human MRPS34 Protein
Catalog #: PEH10768



Summary

Name	RT34_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P82930
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPLFGLGRLV TRKSWLWQHD EPCYWRLTRV RPDYTAQNLD HGKAWGILTF KGKTESEARE IEHVMYHDWR LVPKHEEEAF TAFTPAPEDS LASVPYPPLL RAMIIAERQK NGDTSTEEP LNVQRIRMEP WDYPAKQEDK

image



Background

Alternative Names	28S ribosomal protein S34, mitochondrial (MRP-S34) (S34mt) (Mitochondrial small ribosomal subunit protein mS34)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use