



Product Name: Recombinant Human SERF2 Protein
Catalog #: PEH10772



Summary

Name	SERF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P84101
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTRGNQRELA RQKNMKKQSD SVKGKRRDDG LSAAARKQRD SEIMQQKQKK ANEKKEEPK

image



Background

Alternative Names	Small EDRK-rich factor 2 (Gastric cancer-related protein VRG107) (Protein 4F5-related) (4F5rel) (h4F5rel)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use