


Product Name: Recombinant Human MRPL30 Protein
Catalog #: PEH10786



Summary

Name	RM30_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TCC3
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IPEKVFQASP EDHEKYGGDP QNPHKLHIVT RIKSTRRRPY WEKDIKMLG LEKAHTPQVH KNIPSVNAKL KVVKHLIRIK PLKLPQGLPA

image



Background

Alternative Names	39S ribosomal protein L30, mitochondrial (L30mt) (MRP-L30) (39S ribosomal protein L28, mitochondrial) (L28mt) (MRP-L28) (Mitochondrial large ribosomal subunit protein uL30m)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use