


Product Name: Recombinant Human C2orf49 Protein
Catalog #: PEH10792



Summary

Name	ASHWN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BVC5
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TDSELLHPE LLSQEFLLLT LEQKNIAVET DVRV NKDSL T DLYVQHAIPL PQRDL PKNRW GKMM EKKREQ HEIKN ETKRS STVDGLRKRP LIVFDGSSTS TSIKVKKTEN GDNDRLKPPP QASFTSNAFR KLSNSSSSVS PLILSSNLPV NNKTEHNNND AKQN HDLTHR KSPSGPVKSP PLSPVGTTPV KLKRAAPKEE AEAMNNLKPP QAKRKIQHVT

image



Background

Alternative Names Ashwin

Note

Product Name: Recombinant Human C2orf49 Protein
Catalog #: PEH10792



For Research Use Only , Not for Diagnostic Use