


Product Name: Recombinant Human GOLGA8K Protein
Catalog #: PEH10796



Summary

Name	GOG8K_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	D6RF30
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QHS LQRKGEL ESVLSNVMAT QKKKANQLSS RSKARTEWKL EQSMREEALL KVQLTQLKES FQQVQLERDE YSEHLKGERA RWQQRMRKMS QEICTLKKEK QQDMRRVEKL ERSLSKLKNQ MAEPLPPEPP AVPSEVELQH LRKELERVAG ELQAQVKKNQ RISLLNQRQE ERIQEQEERL RKQEERIQEQ HKSLQQLAKP QSVFEPPNNE NKNALQLEQQ VKELQEKLGE EHLEAASQQN QQLTAQLSLM ALPGEGHGEH LDSEGEEAPQ PMPSVPEDLE

image



Background

Alternative Names Golgin subfamily A member 8K

Product Name: Recombinant Human GOLGA8K Protein
Catalog #: PEH10796



Note

For Research Use Only , Not for Diagnostic Use