


Product Name: Recombinant Human GOLGA6D Protein
Catalog #: PEH10798



Summary

Name	GOG6D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0CG33
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQEEMLRQE AQRVREQERL CEQNERLREQ QKTLQEQGER LRKQEQRLRK QEERLRKEEE RLQKQEKRLW DQEERLWKKE ERLQKQEERL ALSQNHKLDK QLAEPQCSFE DLNNEKKSAL QLEQQVKELQ EKLDEEHLEA ASQRNQQLET QLSLVALPGE GDGGQHLDSE EEEAPRPTPN IPEDLESREA

image



Background

Alternative Names Golgin subfamily A member 6D

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GOLGA6D Protein
Catalog #: PEH10798

