


Product Name: Recombinant Human CT45A5 Protein
Catalog #: PEH10799



Summary

Name	CT455_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P0DMU8
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DPETVFKRPR ECDSPSYQKR QRMALLARKQ GAGDSLIAGS AMSKEKKLMT GHAIPPSQLD SQIDDFTGFS KDGMMQKPGS NAPVGGNVTS SFSGDDLECR ETASSPKSQR EINADIKRKL VKELRCVGQK

image



Background

Alternative Names	Cancer/testis antigen family 45 member A5 (Cancer/testis antigen 45-5) (Cancer/testis antigen 45A5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use