


Product Name: Recombinant Human CFAP210 Protein
Catalog #: PEH10800



Summary

Name	CC173_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q0VFZ6
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QIKEHEEEEE RRKYEEKDA EEIKRQNALY EIEMRKKLEK KREEMHESRR RFLEHMQDKH IIKAVEQQQQ EEEDEKMRKF IKAKKRLIQM GKEKEAETHR LMEKRRERIH NFLSELLKEK LDNEDMIAR DIAEAEAEWE KREREKDEKN KAELKTIAEY RAIVMKNKEE EERQRKIEAK EQLLAVMKAD QIFWEHEKEK KCKADKEHQE

image



Background

Alternative Names Coiled-coil domain-containing protein 173

Note

Product Name: Recombinant Human CFAP210 Protein
Catalog #: PEH10800



For Research Use Only , Not for Diagnostic Use