


Product Name: Recombinant Human TXNRD3NB Protein
Catalog #: PEH10802



Summary

Name	TR3N_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6F5E7
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPELKAEQQM PLEPVRARLS VGLACCCSHT TAEASSLEHG DKVFGQGFPS PLEEIKRLLK ISRALQARSV PSTQEKAKCL SGEPGQPEGK GQETYPGPGK VEGKAEPAMR

image



Background

Alternative Names	Protein TXNRD3NB (Thioredoxin reductase 2 intronic transcript 1) (Thioredoxin reductase 3 intronic transcript 1) (Thioredoxin reductase 3 neighbor gene protein) (TXNRD3 neighbor gene protein) (Thioredoxin reductase 3 new transcript 1)
--------------------------	--

Note



Product Name: Recombinant Human TXNRD3NB Protein
Catalog #: PEH10802

For Research Use Only , Not for Diagnostic Use