


**Product Name: Recombinant Human OGFOD2 Protein**  
**Catalog #: PEH10803**



## Summary

<b>Name</b>	OGFD2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6N063
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RFRGEQQLRR DYGPILRSRG CVSAKDFQQL LAELEQEVER RQRLGQESAA RKALIASYH PARPEVYDSL QDAALAPEFL AVTEYSVSPD ADLKGLLQRL ETVSEEKRIY RVPVFTAPFC QALLELEHF EQSDMPKGRP NTMNNYGVL HELGLDEPLM TPLRERFLQP LMALLYPCG GGRLDSHRAF VVKYAPGQDL ELGCHYDNAE LTLNVALGKV FTGGALYFGG LFAQPTALTE PLEVEHVVGQ GVLHRGGQLH GARPLGTGER WNLVVWLRAS AVRNSLCPMC CREPDLVDDE GFGDGFTRREE

## image



## Background

**Alternative Names** 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 2 (EC 1.14.11.-)

**Product Name: Recombinant Human OGFOD2 Protein**  
**Catalog #: PEH10803**

---



**Note**

For Research Use Only , Not for Diagnostic Use