

Product Name: Recombinant Human CXorf38 Protein
Catalog #: PEH10807



Summary

Name	CX038_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TB03
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKVSSTWLRD FQMKIQNFLN EFKNIP EIVA VYSRIEQLLT SDWAVHIPEE DQRDGCECEM GTYLSSEQVN EIEMQLLKEK LQEIYLQAEQ QEVLPEELSN RLEVVK EFLR NNEDLRNGLT EDMQKLD SLC LHQKLD SQEP

image



Background

Alternative Names Uncharacterized protein CXorf38

Note

For Research Use Only , Not for Diagnostic Use