


Product Name: Recombinant Human METTL26 Protein
Catalog #: PEH10809



Summary

Name	MTL26_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96S19
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MLVAAA AERN KDPILHVL RQ YLDP AQRGV R VLEV ASGSGQ HAAH FARAFP LAEW QPSD VD QRCL DSIA AT TQAQ GLTN VK APLH LDVT WG WEHW GGIL PQ SLD LLLC INM AHVS PLR CTE GLFRA AGH LL KPRALL ITYG PYA INGKISP QSNV DFDL ML RCRN PEWGL R DTALL EDLGK ASG LLLERM V DMPANN KCLI FRKN

image



Background

Alternative Names Methyltransferase-like 26

Note

Product Name: Recombinant Human METTL26 Protein
Catalog #: PEH10809



For Research Use Only , Not for Diagnostic Use